

Anatomage

System Requirements

Updated: November 2017

Summary

IMPORTANT: InvivoMac is not compatible with i-Phone, i-Pad, i-Pad Mini, and Apple Watch. InvivoMac may not be compatible with systems older than 2011.

	Minimum	Recommended
System	2011 editions of <ul style="list-style-type: none">MacbookiMacMac Mini <p><i>with Intel HD 5000 or higher</i></p>	<ul style="list-style-type: none">Latest editions of the Macbook Pro, Macbook Air, Macbook, iMac, or Mac Pro

Detailed Requirements

	Minimum	Recommended
CPU	Intel Core i5	Intel Core i7 or comparable multi-core processor
Graphics Card	Intel HD 5000 <i>or</i> ATI Radeon HD 4650 <i>or</i> Nvidia GeForce 9800 GT Note: Lower-tier graphics cards may run the software but exhibit lower performance.	ATI Radeon HD 6750 <i>or</i> Intel Iris Plus 640 <i>or</i> comparable
RAM	4GB	8GB
Hard Disk	100 GB	500 GB
OS	OS X v 10.10 Yosemite	

Note: Please check the system preferences for the mouse or trackpad. Right-click operations in the PC version are performed by the Mac's secondary click input. Devices such as the Apple Mouse or the Trackpad may impede software navigation. Keyboard commands may also vary depending on the keyboard manufacturer. Some software interfaces may not be displayed correctly if the resolution is set lower than 1280x800.

For more information, please contact Anatomage Tech Support (408) 885-1474.